

SPANDA-SLEEVESTM

Fitting Guidelines

Spanda-Sleeves[™] are intended to fit as a comfort product. Actual fit will depend on individual anthropometric measurements and patient comfort level. Spanda-Sleeves will stretch up to 100% in size (combination of both width and length).

Product	Length	Fits Limb Circumference
# MTPS26000S	15-1/2"	8-1/2"
# MTPS26000	16-1/2"	11"
# MTPS26000L	18"	11"
# MTPS26000XL	20"	13"
# MTPS26001	18-1/4"	11"

CAUTION! The higher the degree of stretching, the tighter the fit and the higher amount of compression on the patient's skin and muscle.

ARM APPLICATION INSTRUCTIONS:

- 1. Beginning with upper arm opening, roll inside edge of Spanda-Sleeve outwards; continue to roll until you reach the thumb opening.
- 2. Slip the arm opening of the Spanda-Sleeve over the patient's hand, fitting the patient's thumb through the thumb hole.
- 3. Unroll the Spanda-Sleeve back up over the arm. Adjust to desired placement.
- 4. If desired, the patient's thumb can be removed from the thumb opening and the Spanda-Sleeve "rolled up" on the wrist to uncover the palm.

LEG APPLICATION INSTRUCTIONS:

- 1. Beginning with upper leg opening, roll inside edge of Spanda-Sleeve outwards until you get within 3" of the toe opening.
- 2. Slip the upper leg opening of the Spanda-Sleeve over the patient's foot, pulling it up until the heel "crook" of the Spanda-Sleeve fits snugly over the patient's heel.
- 3. Unroll the Spanda-Sleeve up the leg. Adjust to desired placement.

WARNING: Inspect skin frequently for tears, lacerations, avulsions, abrasions, infection or other damage and determine appropriate treatment. Change product with new or laundered product in accordance with facility policy.