

MT PROTECTOR-SLEEVES™

Your Best Defense for Protection of At-Risk Skin!

LATEX FREE



PSARM1 - Arm Sleeve



PSLEG1 - Leg Sleeve

PROTECTS ARMS & LEGS WITH SOFT DURABLE FABRIC.

Available in White only.

Fabric consists of 93% Cotton & 7% Spandex.

- **Protects fragile skin**
- **Promotes warmth**
- **Stays in place**
- **Easy to Apply over Arms or Legs**

Seamless, knitted tubular sleeves with finished hems on both ends.

Designed to protect fragile skin against tears, scrapes, abrasions and bruises.

Can also be used to protect IV and wound sites.

PRODUCT INFORMATION:

#PSARM1

**Arm Sleeve - Medium Size - 18" length
with thumb cut-out - 25 Pair/Case**

#PSLEG1

**Leg Sleeve - Universal Size - 14" length
25 Pair/Case**

