

Made in  
the USA

# SPANDAGE™

TUBULAR ELASTIC RETAINER NET

**LATEX FREE**



Spandage™ is manufactured solely  
in the U.S.A. US Patent No. 3747374

- Replaces adhesive tape & wrap gauze
- Holds both wet & dry dressings securely
- 11 Cut-to-fit sizes
- A lucite applicator is included in each box of S01 & S02 at no additional cost
- No allergic reactions. No surprises
- Eliminates skin stripping & excoriation caused by tape removal
- No "wet tape" or tape residue removal
- Equal pressure retains dressing without slippage or excessive constriction
- Dressing changes are easier, faster, & more comfortable - less time for staff, less stressful for patients

HCPC'S CODE A6457

Product #	Description /Application
S01	Small finger, toe w/applicator
S02	Adult finger, toe w/applicator
S03	Average hand, multiple fingers, lower arm, wrist, ankle
S04	Average hand, foot, upper arm, elbow
S05	Average/large arm, small leg, knee amputations
S06	Average/large leg, amputations, intravenous arm board
S07	Average size head, cranium, ears, face, neck
S08	Large leg, cranium, ears, face, neck
S09	Average chest, abdomen, breast, shoulder
S10	X-Large chest, abdomen, breast, shoulder
S11	XX-Large chest, abdomen, breast, shoulder

**25 Yards Per Box - 10 Boxes Per Case**

**M**  
**MEDI-TECH**  
International Corporation

(800) 333-0109

www.medi-techintl.com