


SPANDANETTE™

NON-STERILE STOCKINETTE

NOT MADE WITH Natural Rubber Latex

SpandaNette™ provides light padding and skin protection when used under casts or as protection over wound dressings.



PRODUCT INFORMATION

(HCPCS Code A6457)

Non-Sterile Single Ply Cotton

ORDER#	SIZE & COLOR	DESCRIPTION	CASE QTY.
MTSN13225	2" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	12 BXS
MTSN13325	3" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	12 BXS
MTSN13425	4" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	10 BXS
MTSN13525	5" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	6 BXS
MTSN13625	6" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	6 BXS
MTSN13825	8" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	4 BXS
MTSN131025	10" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	4 BXS
MTSN131225	12" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	2 BXS
MTSN131425	14" - Natural	Non-Sterile Single Ply Cotton Stockinette 25 Yards	2 BXS

Non-Sterile Bias Cut Cotton

ORDER#	SIZE & COLOR	DESCRIPTION	CASE QTY.
MTSNB14250	2" - Natural	Non-Sterile Bias Cut Stockinette 50 Yards	12 BXS
MTSNB14350	3" - Natural	Non-Sterile Bias Cut Stockinette 50 Yards	12 BXS
MTSNB14450	4" - Natural	Non-Sterile Bias Cut Stockinette 50 Yards	10 BXS
MTSNB14550	5" - Natural	Non-Sterile Bias Cut Stockinette 50 Yards	6 BXS
MTSNB14650	6" - Natural	Non-Sterile Bias Cut Stockinette 50 Yards	6 BXS
MTSNB14850	8" - Natural	Non-Sterile Bias Cut Stockinette 50 Yards	4 BXS

Non-Sterile Synthetic (Polyester) Single Ply

ORDER#	SIZE & COLOR	DESCRIPTION	CASE QTY.
MTSNP15225	2" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	12 BXS
MTSNP15325	3" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	12 BXS
MTSNP15425	4" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	10 BXS
MTSNP15625	6" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	6 BXS
MTSNP15825	8" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	4 BXS
MTSNP151025	10" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	4 BXS
MTSNP151225	12" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	4 BXS
MTSNP151425	14" - White	Non-Sterile Single Ply Polyester Stockinette 25 Yards	2 BXS