


SECURING DRESSINGS IS EASY WITH "SPANDAGE" TUBULAR RETAINER NET!



- NO tape or wrap gauze needed!
- Holds both wet & dry dressings firmly
- Equal pressure retains dressings without slippage or excessive constriction



- Elastic Retainer Net fits all body sites; fingers, toes, arms, legs, head, stomach, chest...



**MT Spandage is Latex-Free
and Made in the USA!**

MT SPANDAGE™


TUBULAR ELASTIC RETAINER NET

PEDIATRIC SIZING CHART

LATEX FREE

NEONATAL / X-SMALL CHILD

MT03	Neonatal Head
MT02	X-Small Head
MT02	X-Small Shoulder
MT05	Neonatal Chest, Back, Perineum
MT04	X-Small Chest, Back, Perineum
MT02	Neonatal Arm & Hand
MT01	X-Small Arm & Hand
MTF1	X-Small Fingers
MT03	Neonatal Thigh
MT02	X-Small Thigh
MT03	Neonatal Leg, Foot
MT02	X-Small Leg Foot
MTF1	X-Small Toes



Made in the USA!

SMALL / MEDIUM CHILD

MT04	Small Head
MT05	Medium Head
MT04	Small Shoulder
MT05	Medium Shoulder
MT06	Small Chest, Back, Axilla, Perineum
MT07	Medium Chest, Back, Axilla, Perineum
MT03	Small Arm & Hand
MT04	Medium Arm & Hand
MT02	Small Hand
MTF1	Small Fingers
MTF2	Medium Fingers
MT05	Small Thigh
MT06	Medium Thigh
MT03	Small Leg, Foot
MT04	Medium Leg, Foot
MTF1	Small Toes
MTF2	Medium Toes

Latex Free Tubular Elastic Dressing Retainer Net contains no rubber to cause latex reactions and is available in sizes which hold dressings securely to a variety of body parts and shapes. The information shown above is provided as an example of size/application.

Always move up in size to hold the dressing in place without constriction.



MEDI-TECH
International Corporation

(800) 333-0109

www.medi-techintl.com